



## PacBio pb-16S-nf pipeline

[VIB - Nucleomics Core, nucleomics@vib.be]


*Fri Jun 07, 2024 - version 1.0*

### General Information

Analyse demultiplexed data obtained from one or more Sequel-IIe 16S amplicon sequencing runs and return user-readable summary data and data objects that the end-user can evaluate further using the QIIME2 toolbox.

The analysis was performed using default pipeline parameters except for the Rarefaction which was set to ‘auto’ in order to better compare separate run experiments (this value is generally optimal for a regular run where most of the samples return more than 8000-10000 HiFi reads). Please refer to the separate **barcode\_QC\_v11.html** document to find out about your sample depth and see which samples were below that arbitrary limit.

**Note:** Some of the default options might not be ideal for every biological question and project. If a new analysis needs be done with other parameters, the user can contact us for a quote or use the pb-16S-nf tool on their own bioinformatics infrastructure <sup>1</sup>.



(<https://github.com/PacificBiosciences/HiFi-16S-workflow>)



more at <http://www.nucleomics.be>

<sup>1</sup><https://github.com/PacificBiosciences/HiFi-16S-workflow>

## Analysis settings

The version of the **pb-16S-nf** pipeline used here was **v0.6** <sup>2</sup>.

**Table 1:** list of analysis arguments accepted by the pb-16S-nf pipeline

Other important options:	SET
-front_p Forward primer sequence. Default to F27. (default: AGRGTTYGATYMTGGCTCAG)	AGRGTTYGATYMTGGCTCAG
-adapter_p Reverse primer sequence. Default to R1492. (default: AAGTCGTAACAAGGTARCY)	AAGTCGTAACAAGGTARCY
-filterQ Filter input reads above this Q value (default: 20).	20
-downsample Limit reads to a maximum of N reads if there are more than N reads (default: off)	OFF
-max_ee DADA2 max_EE parameter. Reads with number of expected errors higher than this value will be discarded (default: 2)	2
-minQ DADA2 minQ parameter. Reads with any base lower than this score will be removed (default: 0)	0
-min_len Minimum length of sequences to keep (default: 1000)	1000
-max_len Maximum length of sequences to keep (default: 1600)	1600
-pooling_method QIIME 2 pooling method for DADA2 denoise see QIIME 2 documentation for more details (default: "pseudo", alternative: "independent")	pseudo
-maxreject max-reject parameter for VSEARCH taxonomy classification method in QIIME 2 (default: 100)	100
-maxaccept max-accept parameter for VSEARCH taxonomy classification method in QIIME 2 (default: 100)	100
-min_asv_totalfreq Total frequency of any ASV must be above this threshold across all samples to be retained. Set this to 0 to disable filtering (default 5)	5
-min_asv_sample ASV must exist in at least min_asv_sample to be retained. Set this to 0 to disable. (default 1)	1
-vsearch_identity Minimum identity to be considered as hit (default 0.97)	0.97
-rarefaction_depth Rarefaction curve "max-depth" parameter. By default the pipeline automatically select a cut-off above the minimum of the denoised reads for >80% of the samples. This cut-off is stored in a file called "rarefaction_depth_suggested.txt" file in the results folder (default: null)	'auto' (*)
-skip_primer_trim Skip all primers trimming (switch off cutadapt and DADA2 primers removal) (default: trim with cutadapt)	trim
-skip_nb Skip Naive-Bayes classification (only uses VSEARCH) (default: false)	
-colorby Columns in metadata TSV file to use for coloring the MDS plot in HTML report (default: condition)	condition=run#
-run_picrust2 Run PICRUST2 pipeline. Note that pathway inference with 16S using PICRUST2 has not been tested systematically (default: false)	FALSE

NOTE: We used mostly default values unless noted by (\*). The list of applied parameters can be found in **nextflow\_reports/parameters.txt**.

## Zymo Control sample

When performing the amplification of the 16S V1V9 amplicon for the customer, we usually include a negative control (buffer) and a positive control from the Zymo mock community <sup>3</sup> (ref:**D6305**). These samples are labelled 'Neg\_ctrl' and 'Pos\_ctrl' in our results and will correspond to different barcode pairs in each experiment (amplified alongside with the customer samples).

A PDF file found on the Zymo site describes the mock sample in details <sup>4</sup>

last edits: *FriJun07, 2024*

<sup>2</sup><https://github.com/PacificBiosciences/HiFi-16S-workflow>

<sup>3</sup><https://zymoresearch.eu/collections/zymbiomics-microbial-community-standards/products/zymbiomics-microbial-community-dna-standard>

<sup>4</sup>[https://files.zymoresearch.com/datasheets/ds1706\\_zymbiomics\\_microbial\\_community\\_standards\\_data\\_sheet.pdf](https://files.zymoresearch.com/datasheets/ds1706_zymbiomics_microbial_community_standards_data_sheet.pdf)